



Prompt gamma imaging of a clinical proton beam using Cd Zn Te Compton cameras with a maximum likelihood expectation maximization technique

Josiah De Klerk

University of Cape Town


Department of Physics

Farshad Safavi, Jerimy Polf, Lei Ren, Sina Mossahebi, Steve Peterson



DEPARTMENT OF PHYSICS
UNIVERSITY OF CAPE TOWN
IYUNIVESITHI YASEKAPA • UNIVERSITEIT VAN KAAPSTAD

M e A S U R e
Metrological and Applied Sciences University Research Unit



PGI
CAPE TOWN

UCT
PROTON
THERAPY
INITIATIVE

Cancer – A treatable problem



Globally

2020

A leading cause of death with 10 million deaths 1 in 6 deaths



World Health Organization



Africa

2022:

1,185,216 cancer cases with 763,843 deaths

100,000 new childhood cancer cases are being diagnosed annually in Africa.



Treatment

Surgery

Immunotherapy

Chemotherapy


Radiation Treatment

Radiation treatment: Proton therapy



Goal: To **deliver** a therapeutic dose to the **tumor** while **minimizing radiation** exposure to surrounding **healthy** tissues



How: Precise treatment planning and accurate dose delivery.



The Far Side – Gary Larson



Radiation treatment: Proton therapy

Goal: To deliver a therapeutic dose to the tumor while minimizing radiation exposure to surrounding healthy tissues



How: Precise treatment planning and accurate dose delivery.



Benefits


The low entry point and sharp distal fall-off i.e., protons stop (Finite range)

The low-dose plateau with low energy transfer after the Bragg peak

Limitations


Beam range uncertainty: tissue density fluctuation, range straggling, large penumbra

Patient and organ movement





Prompt Gamma Imaging (PGI)

A solution to proton range verification




During proton therapy, accelerated protons excite atomic nuclei, which promptly de-excite and emit prompt gamma rays.



**200 MeV
protons
on water
target**

Prompt Gamma Imaging (PGI)



Compton Cameras



Each event in the Compton camera produces a cone-of-origin.

$$\cos \theta_1 = 1 m_e c^2 \left(\frac{1}{E_0} - \frac{1}{E_0 - E_1} \right)$$



The cones are overlapped to produce an image of the source location.




**^{137}Cs point source
using CORE
algorithm from the
UCT Polaris
detector system**

CdZnTe CC Detector Systems

UCT



UMMC Nozzle-mounted



- Four $20 \times 20 \times 10$ mm³ cadmium zinc telluride crystals
- Each crystal arranged in a 2×2 array, pixelated in a 11×11 x- and y- direction and a z- direction depth of interaction

3D Gamma detection

Room Temperature


Good Energy Resolution (0.8%

FWHM at 662 keV)

Compact

Gamma Attenuation

In a nozzle-mounted detector system, reconstruction along the direction of the beam is challenging due to attenuation of the gammas as they travel through the patient.



**Raw
Reconstruction**

**Corrected
Reconstruction**

**Monte Carlo
Simulations of 200
MeV protons incident
on water target of the
Maryland Clinical
System**

Maryland Proton Treatment Centre in Baltimore

1

Feasibility — Can the detector system reconstruct proton depth deposition?

2

Investigate gamma-ray attenuation and range shift.

3

Using the UCT detector system to support the UMMC detector system

150 MeV single-spot beam
Gantry 270°
80 nA · 60 kMU











Image Recon. Pipeline




UCT Detector Results




Total Events	4.99×10^6
Preprocessing	3.79×10^5
CORE Filtering	2.84×10^5
Used	2.84×10^5

Parameter	Phantom	CORE
Dims. (mm)	1280^3	800^3
Grid	64^3	100^3
Voxel (mm)	20^3	8^3


Experimental Setup



UCT Detector Results




Expected Location




UCT Detector Results

Radiation Therapy Plan and Dose






CT

DICOM Image with Overlay Plot: $z = 0.0$ mm




UCT Detector Results

Expected Location : (-148.0, 4.0, 356.0) mm



UMMC Nozzle-mounted Results





Total Events	3.79×10^7
Preprocessing	3.26×10^6
CORE Filtering	3.22×10^6
Used	2.00×10^6




Parameter	Phantom	CORE
Dims. (mm)	1280^3	1200^3
Grid	64^3	150^3
Voxel (mm)	20^3	8^3

Experimental Setup UMMC Nozzle-mounted Results

Experimental Setup




Expected Location




UMMC Nozzle-mounted Results

Radiation Therapy Plan and Dose






CT





UMMC Nozzle-mounted Results


Expected Location : (-92.0, 4.0, 12.0) mm



UMMC Nozzle-mounted Results



UMMC Nozzle-mounted Results




Axis	Max. (mm)	Fit (mm)	FWHM (mm)	σ (mm)
X	-92.00	—	—	—
Y	4.00	-3.98	239.12	101.54
Z	12.00	10.34	274.88	116.73

Distance (Estimated - Actual)


92.87 mm

UMMC 10 mm Range Shift Results

Expected Location : (-108.0, -12.0, 4.0) mm



UMMC 10 mm Range Shift Results




Axis	Max. (mm)	Fit (mm)	FWHM (mm)	σ (mm)
X	-108.00	—	—	—
Y	-12.00	-3.98	236.78	100.55
Z	4.00	11.94	263.63	111.95

Distance (Estimated - Actual)

108.74 mm

Range Shift Comparisons



Metric	No Shift	10 mm Range Shift	Change
Distance Error (mm)	92.87	108.74	+15.87 (+17.1%)
Y FWHM (mm)	239.12	236.78	-2.34 (-1.0%)
Z FWHM (mm)	274.88	263.63	-11.25 (-4.1%)

Conclusion

1

The Quad CdZnTe PGI detector system reconstructed proton depth distributions within a phantom.

2

Gamma attenuation in the nozzle-mounted detector geometry. Sensitivity to beam range changes through a 10 mm proton range shift demonstrated

3

The UCT measurands may require further optimization before it can effectively support the UMMC measurements.

Future work

1

Optimize preprocessing, event filtering and image reconstruction algorithms.

2

Investigate CT-based attenuation correction to compensate for prompt gamma attenuation.

3

Investigate machine learning-based image reconstruction and attenuation correction.

Thank You !